


**Connector type**




**Lead Position Marked**



Standard Snap



type AA Grabber



type RL Grabber



**Type Disignation**

E □ □ □ □ □ □ □ □

5A: 5Lead: A, Standard AHA

3H: 3Lead: H, Standard AHA

3A: 3Lead: A, Standard AHA

5I: 5Lead: I, Standard IEC

3IH: 3Lead: IH, Standard IEC

3I: 3Lead: I, Standard IEC













S: Snap (Patient end termination)

C: Clip (Patient end termination)









Brand code: as 041:Bionet, 040:Biosys









Yoke type A: One-piece, B: Lifelinc, C: Din style, D: AA style, E: VS style











E: ECG Cable with Leadwires











Compatible Brand	Compatible Models No.	INSTRUMENT CONNECTOR	REF NO.	STANDARD	LEADS	PATIENT END TERMINAL	Material	Defibrillation withstand	PRICE CODE	
<b>AIR SHIELDS</b>	<b>(AIR)</b> :HRRM-71-1/2,Porta-Fib I,II,III, System5(new),System6(Neonatal)		EA002S3A	AHA,	3LD:A,	SNAP,	TPU.		B5	
			EA002C3A	AHA,	3LD:A,	CLIP,	TPU.		B6	
<b>ALT</b>	<b>(ALT)</b> :800,800Plus,CX ATL Ultrasound CX200,CX400,Apogee,HDI 9/2000 Interspec Cardiology XL Mark IV,Mark VII ,Ultramark 4/5/7/9		EA002S3I	IEC,	3LD:I,	SNAP,	TPU.		B5	
			EA002C3I	IEC,	3LD:I,	CLIP,	TPU.		B6	
			EA002S5A	AHA,	5LD:A,	SNAP,	TPU.		B6	
			EA002C5A	AHA,	5LD:A,	CLIP,	TPU.		B7	
<b>BCI</b>	<b>(BCI)</b> :3101,3404,6100,6200, 9100,9200,Advisor ,Autocorr		EA002S5I	IEC,	5LD:I,	SNAP,	TPU.		B8	
			EA002C5I	IEC,	5LD:I,	CLIP,	TPU.		B9	
<b>COLIN</b>	<b>(COLIN)</b> :BP308, BP408, CBM3000		EA003S3A	AHA,	3LD:A,	SNAP,	TPU,	1K Ω	B8	
<b>CSI</b>	<b>(CSI)</b> :504US,602-4USP,Poet2,Scholar		EA003C3A	AHA,	3LD:A,	CLIP,	TPU,	1K Ω	B9	
<b>GE MEDICAL</b>	<b>(GMS)</b> :Corometrics 556 Critikon 8100T/Dinamap 8700/9700		EA003S3I	IEC,	3LD:I,	SNAP,	TPU,	1K Ω	B8	
<b>IVY</b>	<b>(IVY)</b> :101 (before 1989),402, 700		EA003C3I	IEC,	3LD:I,	CLIP,	TPU,	1K Ω	B9	
<b>LSI</b>	<b>(LSI)</b> :LS5 , 261		EA003S5A	AHA,	5LD:A,	SNAP,	TPU,	1K Ω	B9	
<b>MDE</b>	<b>(MDE)</b> :Escort 100, 101,102,200, 300 Prism, 360SL, 450SL,521A PortaPak 80,PortaPak 90/1000		EA003C5A	AHA,	5LD:A,	CLIP,	TPU,	1K Ω	B0	
			EA003S5I	IEC,	5LD:I,	SNAP,	TPU,	1K Ω	B9	
				EA003C5I	IEC,	5LD:I,	CLIP,	TPU,	1K Ω	B0
<b>PHYSIO-CONTROL (MEDTRONIC)</b>	<b>(MPC)</b> : Lifepak 6 , VSM-1ESF (w/resistors)									
<b>NELLCOR</b>	<b>(NELLCOR)</b> :N-200/250 C-LOCK N-1000,NPB-3900,NPB-4000									
<b>NIHON KOHDEN</b>	<b>(NK)</b> :BSM-8301A/8302A/8500A/8502A OEC-6102A (s/n>20156)/7102A/8108A OEC-6105A (s/n>20226),TEC-7100A TEC-5200A(s/n>20421),TEC-7200A/7300A;		EA020S3A	AHA,	3LD:A,	SNAP,	TPU.		B5	
			EA020C3A	AHA,	3LD:A,	CLIP,	TPU.		B6	
			EA020S3I	IEC,	3LD:I,	SNAP,	TPU.		B5	
			EA020C3I	IEC,	3LD:I,	CLIP,	TPU.		B6	
			EA020S5A	AHA,	5LD:A,	SNAP,	TPU.		B6	
			EA020C5A	AHA,	5LD:A,	CLIP,	TPU.		B7	
			EA020S5I	IEC,	5LD:I,	SNAP,	TPU.		B8	
			EA020C5I	IEC,	5LD:I,	CLIP,	TPU.		B9	
			EA021S3A	AHA,	3LD:A,	SNAP,	TPU,	1K Ω	B8	
			EA021C3A	AHA,	3LD:A,	CLIP,	TPU,	1K Ω	B9	
<b>PACE-TECH</b>	<b>(PACE-TECH)</b> :Minipack 300 Minipack 3000/3100,MiniMax4000CL Vitalmax 800 Plus/2100/2200/4000/4100 Vitalmax/4000CL/4000 Modular/4100CL		EA021S3I	IEC,	3LD:I,	SNAP,	TPU,	1K Ω	B8	
			EA021C3I	IEC,	3LD:I,	CLIP,	TPU,	1K Ω	B9	
<b>PHILIPS</b>	<b>(PMS)</b> :M3921A A1,M3923A A1 M3925AA1,M3927AA3,M3928A3, M3929AA3		EA021S5A	AHA,	5LD:A,	SNAP,	TPU,	1K Ω	B9	
			EA021C5A	AHA,	5LD:A,	CLIP,	TPU,	1K Ω	B0	
<b>PIONEER</b>	<b>(PIONEER)</b> (American Optical):26125 33167, 5260 , 5330 , Pulsar 4		EA021S5I	IEC,	5LD:I,	SNAP,	TPU,	1K Ω	B9	
<b>PPG</b>	<b>(PPG)</b> (Datamedix, Becton-Dickinson) Guardian 1001/N1001		EA021C5I	IEC,	5LD:I,	CLIP,	TPU,	1K Ω	B0	
<b>SPACELABS</b>	<b>(BURDICK)</b> : Medic3/4/5/4045Defib 500/600/900 Series (w/o respiration) 90300/90400/90600/90700 Series TEK-208/408/511/512/514/521 PC Bedside , PCMS									
<b>TEKTRONIX</b>	<b>(Vitatek)</b> :4045Defib/208/400/500Series									
<b>WELCH ALLYN</b>	<b>(Protocol)</b> :Propaq100/200Series,Atlas Propaq Encore, Propaq CS,Ultra Smart									
<b>ZOLL</b>	<b>(ZOLL)</b> :D900, PD1200, PD1400									
<b>GOLDWAY</b>	<b>GOLDWAY</b> : UT4000A, UT4000A, UT4000Apro,UT6000A, UT4000B, UT4000C, UT4000F, UT6000F, UT4000Fpro, UT4800									
<b>Biolight</b>										
<b>Guoteng</b>										

















Compatible Brand	Compatible Models No.	INSTRUMENT CONNECTOR	REF NO.	STANDARD	LEADS	PATIENT END TERMINAL	Material	Defibrillation withstand	PRICE CODE
<b>BIONET KOREA</b>	BM3		EA041S5A	AHA,	5LD:A,	SNAP,	TPU,	1KQ	C1
			EA041C5A	AHA,	5LD:A,	CLIP,	TPU,	1KQ	C2
			EA041S3A	AHA,	3LD:A,	SNAP,	TPU,	1KQ	B9
			EA041C3A	AHA,	3LD:A,	CLIP,	TPU,	1KQ	B0
			EA041S5I	IEC,	5LD:I,	SNAP,	TPU,	1KQ	C1
			EA041C5I	IEC,	5LD:I,	CLIP,	TPU,	1KQ	C2
			EA041S3I	IEC,	3LD:I,	SNAP,	TPU,	1KQ	B9
			EA041C3I	IEC,	3LD:I,	CLIP,	TPU,	1KQ	B0
<b>BIOSYS KOREA</b>	BPM-103		EA040S5A	AHA,	5LD:A,	SNAP,	TPU,	1KQ	C1
			EA040C5A	AHA,	5LD:A,	CLIP,	TPU,	1KQ	C2
			EA040S3A	AHA,	3LD:A,	SNAP,	TPU,	1KQ	B9
			EA040C3A	AHA,	3LD:A,	CLIP,	TPU,	1KQ	B0
			EA040S5I	IEC,	5LD:I,	SNAP,	TPU,	1KQ	C1
			EA040C5I	IEC,	5LD:I,	CLIP,	TPU,	1KQ	C2
			EA040S3I	IEC,	3LD:I,	SNAP,	TPU,	1KQ	B9
			EA040C3I	IEC,	3LD:I,	CLIP,	TPU,	1KQ	B0
<b>BRUKER (ODAM)</b>	SM784, SM785 with 12-pin input connector		EA017S3H	AHA,	3LD:H,	SNAP,	TPU,		C3
			EA017C3H	AHA,	3LD:H,	CLIP,	TPU,		C4
			EA017S3IH	IEC,	3LD:IH,	SNAP,	TPU,		C3
			EA017C3IH	IEC,	3LD:IH,	CLIP,	TPU,		C4
<b>COLIN</b>	BP88 / BP306		EA062S5A	AHA,	5LD:A,	SNAP,	TPU,		D5
			EA062C5A	AHA,	5LD:A,	CLIP,	TPU,		D6
			EA062S3A	AHA,	3LD:A,	SNAP,	TPU,		D3
			EA062C3A	AHA,	3LD:A,	CLIP,	TPU,		D4
			EA062S5I	IEC,	5LD:I,	SNAP,	TPU,		D5
			EA062C5I	IEC,	5LD:I,	CLIP,	TPU,		D6
			EA062S3I	IEC,	3LD:I,	SNAP,	TPU,		D3
			EA062C3I	IEC,	3LD:I,	CLIP,	TPU,		D4
	BP88S		EA061S5A	AHA,	5LD:A,	SNAP,	TPU,		D2
			EA061C5A	AHA,	5LD:A,	CLIP,	TPU,		D3
			EA061S3A	AHA,	3LD:A,	SNAP,	TPU,		D1
			EA061C3A	AHA,	3LD:A,	CLIP,	TPU,		D2
			EA061S5I	IEC,	5LD:I,	SNAP,	TPU,		D2
			EA061C5I	IEC,	5LD:I,	CLIP,	TPU,		D3
			EA061S3I	IEC,	3LD:I,	SNAP,	TPU,		D1
			EA061C3I	IEC,	3LD:I,	CLIP,	TPU,		D2
<b>Creative科瑞康 Spurt世博达 Biolight宝莱特</b>	Creative:PC9000S,PC9000,UP6000, Biolight: M12		EA046S5A	AHA,	5LD:A,	SNAP,	TPU,		C1
			EA046C5A	AHA,	5LD:A,	CLIP,	TPU,		C2
			EA046S3A	AHA,	3LD:A,	SNAP,	TPU,		B9
			EA046C3A	AHA,	3LD:A,	CLIP,	TPU,		B0
			EA046S5I	IEC,	5LD:I,	SNAP,	TPU,		C1
			EA046C5I	IEC,	5LD:I,	CLIP,	TPU,		C2
			EA046S3I	IEC,	3LD:I,	SNAP,	TPU,		B9
			EA046C3I	IEC,	3LD:I,	CLIP,	TPU,		B0
<b>Corpuls</b>	Corpuls 08/16, 0816M 		EA076C3I	IEC,	3LD:IH,	CLIP,	TPU,		E0

Compatible Brand	Compatible Models No.	INSTRUMENT CONNECTOR	REF NO.	STANDARD	LEADS	PATIENT END TERMINAL	Material	Defibrillation withstand	PRICE CODE			
<b>DATASCOPE</b>	2000,2001,2002,2100,2200,3000, 3000A ,Ballon Pump, M/D3 w/lead select,M/D3A,M/D4, Passport,Passport IL,Passport XG, Passport XL,Passport 2, Passport 2LT,Spectrum System 90, System 90T , System 96, System 97		EA002S5A	AHA,	5LD:A,	SNAP,	TPU,		B9			
			EA002C5A	AHA,	5LD:A,	CLIP,	TPU,		B0			
			EA002S3A	AHA,	3LD:A,	SNAP,	TPU,		B7			
			EA002C3A	AHA,	3LD:A,	CLIP,	TPU,		B8			
			EA002S5I	IEC,	5LD:I,	SNAP,	TPU,		B9			
			EA002C5I	IEC,	5LD:I,	CLIP,	TPU,		B0			
			EA002S3I	IEC,	3LD:I,	SNAP,	TPU,		B7			
			EA002C3I	IEC,	3LD:I,	CLIP,	TPU,		B8			
<b>DATEX</b>	AS/3, CS/3, S/5, Cardiocap 1, Cardiocap 2 , Cardiocap 5 , Light Monitor , Satliteplus		EA010S5A	AHA,	5LD:A,	SNAP,	TPU,		C9			
			EA010C5A	AHA,	5LD:A,	CLIP,	TPU,		C0			
			EA010S3A	AHA,	3LD:A,	SNAP,	TPU,		C7			
			EA010C3A	AHA,	3LD:A,	CLIP,	TPU,		C8			
			EA010S5I	IEC,	5LD:I,	SNAP,	TPU,		C9			
			EA010C5I	IEC,	5LD:I,	CLIP,	TPU,		C0			
			EA010S3I	IEC,	3LD:I,	SNAP,	TPU,		C7			
			EA010C3I	IEC,	3LD:I,	CLIP,	TPU,		C8			
<b>DRAGER SIEMENS</b>	Drager UM1/3/3.1, RM1, PM8010, PM8014, Babyguard 8000, Cato, Babylog8000/8000 plus, PB8800 Serie 8000, Cicero, Parameterbox,		EA057S3A	AHA,	3LD:A,	SNAP,	TPU,		D0			
			EA057C3A	AHA,	3LD:A,	CLIP,	TPU,		E1			
			EA057S3I	IEC,	3LD:I,	SNAP,	TPU,		D0			
			EA057C3I	IEC,	3LD:I,	CLIP,	TPU,		E1			
<b>FUKUDA DENSHI</b>	DATASCOPE:Expert FUKUDA DENSHI:DS5100E, DS5300W, DS5400,		EA068S5A	AHA,	5LD:A,	SNAP,	TPU, 1KΩ		D4			
			EA068C5A	AHA,	5LD:A,	CLIP,	TPU, 1KΩ		D5			
			EA068S3A	AHA,	3LD:A,	SNAP,	TPU, 1KΩ		D3			
			EA068C3A	AHA,	3LD:A,	CLIP,	TPU, 1KΩ		D4			
			EA068S5I	IEC,	5LD:I,	SNAP,	TPU, 1KΩ		D4			
			EA068C5I	IEC,	5LD:I,	CLIP,	TPU, 1KΩ		D5			
			EA068S3I	IEC,	3LD:I,	SNAP,	TPU, 1KΩ		D3			
			EA068C3I	IEC,	3LD:I,	CLIP,	TPU, 1KΩ		D4			
<b>PETAS</b>			EA064S3H	AHA,	3LD:H,	SNAP,	TPU,		C8			
			EA064C3H	AHA,	3LD:H,	CLIP,	TPU,		C9			
			EA064S3IH	IEC,	3LD:IH,	SNAP,	TPU,		C8			
			EA064C3IH	IEC,	3LD:IH,	CLIP,	TPU,		C9			
<b>GE-MEDICAL MARQUETTE</b>	GEMS:Dash PRO 4000, Dash PRO 3000, Dash PRO 2000, Dash 1000, Eagle Systems, MAC-Lab, Solar 9500, Solar 8000. Solar 3000, Tram 100, Tram 200, Tram 300, Tram 400, Tram 500, Tram 600, Tram 800, Unity Systems		EA007S5A	AHA,	5LD:A,	SNAP,	TPU,		C7			
			EA007C5A	AHA,	5LD:A,	CLIP,	TPU,		C8			
			EA007S3A	AHA,	3LD:A,	SNAP,	TPU,		C6			
			EA007C3A	AHA,	3LD:A,	CLIP,	TPU,		C7			
			EA007S5I	IEC,	5LD:I,	SNAP,	TPU,		C7			
			EA007C5I	IEC,	5LD:I,	CLIP,	TPU,		C8			
			EA007S3I	IEC,	3LD:I,	SNAP,	TPU,		C6			
			EA007C3I	IEC,	3LD:I,	CLIP,	TPU,		C7			
	PRO1000		EA021S3A	AHA,	3LD:A,	SNAP,	TPU, 1KΩ		B7			
			EA021C3A	AHA,	3LD:A,	CLIP,	TPU, 1KΩ		B8			
			EA021S3I	IEC,	3LD:I,	SNAP,	TPU, 1KΩ		B7			
			EA021C3I	IEC,	3LD:I,	CLIP,	TPU, 1KΩ		B8			
			<b>GE-MEDICAL</b>	GE MEDICAL SYSTEMS:2000, 2020, 2110, 2110N, 2120, 2210, 2230, 2230N, 2310, 2330, 2330N, 7000 Series, Tram AR		EA069S5A	AHA,	5LD:A,	SNAP,	TPU,		C9
						EA069C5A	AHA,	5LD:A,	CLIP,	TPU,		C0
EA069S3A	AHA,	3LD:A,				SNAP,	TPU,		C7			
EA069C3A	AHA,	3LD:A,				CLIP,	TPU,		C8			
EA069S5I	IEC,	5LD:I,				SNAP,	TPU,		C9			
EA069C5I	IEC,	5LD:I,				CLIP,	TPU,		C0			
EA069S3I	IEC,	3LD:I,				SNAP,	TPU,		C7			
EA069C3I	IEC,	3LD:I,				CLIP,	TPU,		C8			

Compatible Brand	Compatible Models No.	INSTRUMENT CONNECTOR	REF NO.	STANDARD	LEADS	PATIENT END TERMINAL	Material	Defibrillation withstand	PRICE CODE				
<b>HELLIGE</b>	Hellige, Servomed 4 , Generation, SCP 851, SCP 852 , SCP 844, SCS 804,, CardioServ , VICOM-sm		EA011S5A	AHA,	5LD:A,	SNAP,	TPU,		D4				
			EA011C5A	AHA,	5LD:A,	CLIP,	TPU,		D5				
			EA011S3A	AHA,	3LD:A,	SNAP,	TPU,		D3				
			EA011C3A	AHA,	3LD:A,	CLIP,	TPU,		D4				
			EA011S5I	IEC,	5LD:I,	SNAP,	TPU,		D4				
			EA011C5I	IEC,	5LD:I,	CLIP,	TPU,		D5				
			EA011S3I	IEC,	3LD:I,	SNAP,	TPU,		D3				
			EA011C3I	IEC,	3LD:I,	CLIP,	TPU,		D4				
<b>KONTRON</b>	7000, 7135, 7135B, 7141, 7210 7250, 7251, K2000, KAAT, Kolormon, Micromon, Minimon, Supermon , Vitalmon		EA031S5A	AHA,	5LD:A,	SNAP,	TPU,		C5				
			EA031C5A	AHA,	5LD:A,	CLIP,	TPU,		C6				
			EA031S3A	AHA,	3LD:A,	SNAP,	TPU,		C3				
			EA031C3A	AHA,	3LD:A,	CLIP,	TPU,		C4				
			EA031S5I	IEC,	5LD:I,	SNAP,	TPU,		C5				
			EA031C5I	IEC,	5LD:I,	CLIP,	TPU,		C6				
			EA031S3I	IEC,	3LD:I,	SNAP,	TPU,		C3				
			EA031C3I	IEC,	3LD:I,	CLIP,	TPU,		C4				
		463		EA032S3A	AHA,	3LD:A,	SNAP,	TPU,		C9			
				EA032C3A	AHA,	3LD:A,	CLIP,	TPU,		C0			
				EA032S3I	IEC,	3LD:I,	SNAP,	TPU,		C9			
				EA032C3I	IEC,	3LD:I,	CLIP,	TPU,		C0			
				<b>L&amp;T (INDIA )</b>	DLanet50		EA030S3H	AHA,	3LD:H,	SNAP,	TPU, 1KΩ		C5
							EA030C3H	AHA,	3LD:H,	CLIP,	TPU, 1KΩ		C6
EA030S3IH	IEC,	3LD:IH,	SNAP,				TPU, 1KΩ		C5				
EA030C3IH	IEC,	3LD:IH,	CLIP,				TPU, 1KΩ		C6				
<b>MEK</b>	MP1000,MP600,MP500		EA029S5A	AHA,	5LD:A,	SNAP,	TPU,		B9				
			EA029C5A	AHA,	5LD:A,	CLIP,	TPU,		B0				
			EA029S3A	AHA,	3LD:A,	SNAP,	TPU,		B7				
			EA029C3A	AHA,	3LD:A,	CLIP,	TPU,		B8				
			EA029S5I	IEC,	5LD:I,	SNAP,	TPU,		B9				
			EA029C5I	IEC,	5LD:I,	CLIP,	TPU,		B0				
			EA029S3I	IEC,	3LD:I,	SNAP,	TPU,		B7				
			EA029C3I	IEC,	3LD:I,	CLIP,	TPU,		B8				
	M160246 M360697	172		EA033S5A	AHA,	5LD:A,	SNAP,	TPU,		C4			
				EA033C5A	AHA,	5LD:A,	CLIP,	TPU,		C5			
				EA033S3A	AHA,	3LD:A,	SNAP,	TPU,		C3			
				EA033C3A	AHA,	3LD:A,	CLIP,	TPU,		C4			
				EA033S5I	IEC,	5LD:I,	SNAP,	TPU,		C4			
				EA033C5I	IEC,	5LD:I,	CLIP,	TPU,		C5			
			EA072S3A	AHA,	3LD:A,	SNAP,	TPU,		C4				
			EA072S5A	AHA,	5LD:A,	SNAP,	TPU,		C3				
<b>MENNEN</b>	Horizon 250,3000,450, Microsolo,936 Series,425/101 Preamp,Cardio-Pak 2000,740, 741,742		EA059S3A	AHA,	3LD:A,	SNAP,	TPU, 1KΩ		C4				
			EA059C3A	AHA,	3LD:A,	CLIP,	TPU, 1KΩ		C5				
			EA059S3I	IEC,	3LD:I,	SNAP,	TPU, 1KΩ		C4				
			EA059C3I	IEC,	3LD:I,	CLIP,	TPU, 1KΩ		C5				
	(MENNEN): 900,1100,1110, Horizon 2000,2100,2200,5200 Horizon 9000XL,MR1200,1270 Horizon 1280,1300,Horizon 1380	175		EA065S5A	AHA,	5LD:A,	SNAP,	TPU,		D9			
				EA065C5A	AHA,	5LD:A,	CLIP,	TPU,		D0			
				EA065S3A	AHA,	3LD:A,	SNAP,	TPU,		D7			
				EA065C3A	AHA,	3LD:A,	CLIP,	TPU,		D8			
				EA065S5I	IEC,	5LD:I,	SNAP,	TPU,		D9			
				EA065C5I	IEC,	5LD:I,	CLIP,	TPU,		D0			
			EA065S3I	IEC,	3LD:I,	SNAP,	TPU,		D7				
			EA065C3I	IEC,	3LD:I,	CLIP,	TPU,		D8				

Compatible Brand	Compatible Models No.	INSTRUMENT CONNECTOR	REF NO.	STANDARD	LEADS	PATIENT END TERMINAL	Material	Defibrillation withstand	PRICE CODE			
<b>MINDRAY</b>	MINDRAY: PM9000,PM8000,PM7000 with 6-pin input connector  <b>MEC1000, MEC1200</b>		EA020S3A	AHA,	3LD:A,	SNAP,	TPU,		B5			
			EA020C3A	AHA,	3LD:A,	CLIP,	TPU,		B6			
			EA020S3I	IEC,	3LD:I,	SNAP,	TPU,		B5			
			EA020C3I	IEC,	3LD:I,	CLIP,	TPU,		B6			
			EA020S5A	AHA,	5LD:A,	SNAP,	TPU,		B6			
			EA020C5A	AHA,	5LD:A,	CLIP,	TPU,		B7			
			EA020S5I	IEC,	5LD:I,	SNAP,	TPU,		B8			
			EA020C5I	IEC,	5LD:I,	CLIP,	TPU,		B9			
			BeneView (T5, T6, T8) iPM series  可兼容原装0010-30-42719五导电缆  可兼容原装0010-30-42720三导电缆		EA050S5A	AHA,	5LD:A,	SNAP,	TPU,		C9	
	EA050C5A	AHA,			5LD:A,	CLIP,	TPU,		C0			
	EA050S3A	AHA,			3LD:A,	SNAP,	TPU,		C7			
	EA050C3A	AHA,			3LD:A,	CLIP,	TPU,		C8			
	EA050S5I	IEC,			5LD:I,	SNAP,	TPU,		C9			
	EA050C5I	IEC,			5LD:I,	CLIP,	TPU,		C0			
	EA050S3I	IEC,			3LD:I,	SNAP,	TPU,		C7			
	EA050C3I	IEC,			3LD:I,	CLIP,	TPU,		C8			
	PM5000/6000 (12孔插座的监护仪)				EA060S5A	AHA,	5LD:A,	SNAP,	TPU,		C3	
			EA060C5A	AHA,	5LD:A,	CLIP,	TPU,		C3			
			EA060S3A	AHA,	3LD:A,	SNAP,	TPU,		C3			
			EA060C3A	AHA,	3LD:A,	CLIP,	TPU,		C3			
			EA060S5I	IEC,	5LD:I,	SNAP,	TPU,		C3			
			EA060C5I	IEC,	5LD:I,	CLIP,	TPU,		C3			
			EA060S3I	IEC,	3LD:I,	SNAP,	TPU,		C3			
			EA060C3I	IEC,	3LD:I,	CLIP,	TPU,		C3			
			<b>NEC</b>	3M12		EA053S5A	AHA,	5LD:A,	SNAP,	TPU,		D4
	EA053C5A	AHA,				5LD:A,	CLIP,	TPU,		D5		
	EA053S3A	AHA,				3LD:A,	SNAP,	TPU,		D3		
EA053C3A	AHA,	3LD:A,				CLIP,	TPU,		D4			
EA053S5I	IEC,	5LD:I,				SNAP,	TPU,		D4			
EA053C5I	IEC,	5LD:I,				CLIP,	TPU,		D5			
EA053S3I	IEC,	3LD:I,				SNAP,	TPU,		D3			
EA053C3I	IEC,	3LD:I,				CLIP,	TPU,		D4			
PB1916E		EA054S3A				AHA,	3LD:A,	SNAP,	TPU,		E9	
		EA054C3A		AHA,	3LD:A,	CLIP,	TPU,		E0			
		EA054S3I		IEC,	3LD:I,	SNAP,	TPU,		E9			
		EA054C3I		IEC,	3LD:I,	CLIP,	TPU,		E0			
		<b>NIHON KOHDEN</b>		OEC-6102A (s/n<20156), OEC-6105A (s/n<20226), OEC-6120A/J/K BSM-7105K		EA037S5A	AHA,	5LD:A,	SNAP,	TPU, 1KΩ		C9
						EA037C5A	AHA,	5LD:A,	CLIP,	TPU, 1KΩ		C0
EA037S3A	AHA,					3LD:A,	SNAP,	TPU, 1KΩ		C7		
EA037C3A	AHA,		3LD:A,			CLIP,	TPU, 1KΩ		C8			
EA037S5I	IEC,		5LD:I,			SNAP,	TPU, 1KΩ		C9			
EA037C5I	IEC,		5LD:I,			CLIP,	TPU, 1KΩ		C0			
EA037S3I	IEC,		3LD:I,			SNAP,	TPU, 1KΩ		C7			
EA037C3I	IEC,		3LD:I,			CLIP,	TPU, 1KΩ		C8			
BSM-7102A,OMP-7201C, OMP-7203C, TEC-5200A (s/n<20421)			EA038S3A			AHA,	3LD:A,	SNAP,	TPU,		D7	
			EA038C3A	AHA,	3LD:A,	CLIP,	TPU,		D8			
			EA038S3I	IEC,	3LD:I,	SNAP,	TPU,		D7			
TEC-7621C			EA075S3A	AHA,	3LD:A,	SNAP,	TPU, 200KΩ		D7			
		EA075C3A	AHA,	3LD:A,	CLIP,	TPU, 200KΩ		D8				
		EA075S3I	IEC,	3LD:I,	SNAP,	TPU, 200KΩ		D7				
		EA075C3I	IEC,	3LD:I,	CLIP,	TPU, 200KΩ		D8				
TEC-7621K, TEC-7631K		EA074S3A	AHA,	3LD:A,	SNAP,	TPU, 100KΩ		D7				
		EA074C3A	AHA,	3LD:A,	CLIP,	TPU, 100KΩ		D8				
		EA074S3I	IEC,	3LD:I,	SNAP,	TPU, 100KΩ		D7				
		EA074C3I	IEC,	3LD:I,	CLIP,	TPU, 100KΩ		D8				
<b>PETAS</b>	KMA160R,165R,260R,K-250,300...		EA063S3H	AHA,	3LD:H,	SNAP,	TPU,		B9			
			EA063C3H	AHA,	3LD:H,	CLIP,	TPU,		B0			
			EA063S3IH	IEC,	3LD:IH,	SNAP,	TPU,		B9			
			EA063C3IH	IEC,	3LD:IH,	CLIP,	TPU,		B0			

Compatible Brand	Compatible Models No.	INSTRUMENT CONNECTOR	REF NO.	STANDARD	LEADS	PATIENT END TERMINAL	Material	Defibrillation withstand	PRICE CODE			
<b>PHILIPS</b>	43100A,43110A,43120A,78352C, 78354C,78834C,862474 C3, M1001A,M1002A,M1165A, M1166A,M1167A,M1175A,M1177A ,M1205A,M1264A,M1275A Transport, M1280A,M1722A/B,M1723A/B, M2475B,M2601A,M3000A,M3001A, M3046A M2/M3/M4,M4735A	018 	EA023S5A	AHA,	5LD:A,	SNAP,	TPU,	1KΩ	C9			
			EA023C5A	AHA,	5LD:A,	CLIP,	TPU,	1KΩ	C0			
			EA023S3A	AHA,	3LD:A,	SNAP,	TPU,	1KΩ	C7			
			EA023C3A	AHA,	3LD:A,	CLIP,	TPU,	1KΩ	C8			
			EA023S5I	IEC,	5LD:I,	SNAP,	TPU,	1KΩ	C9			
			EA023C5I	IEC,	5LD:I,	CLIP,	TPU,	1KΩ	C0			
			EA023S3I	IEC,	3LD:I,	SNAP,	TPU,	1KΩ	C7			
			EA023C3I	IEC,	3LD:I,	CLIP,	TPU,	1KΩ	C8			
	43100A/43110A/43110M/43120A/MC Defib, 43200A Portable, M3500B,77000AC/77020A,77010AC Sonos/AC/AR,77025A,77030A,77035A ,77040A,78208A,78210A,78330A,78331A,78332A,78333A,78341A,78342A,78345A,78346A,78351T,78352A/C,78353A/B,78354A,78532A/B,78534A/B/C, 78619A,78620A,78660A/B,78661A/B,78670A,78671A,M1722A/B,M1723A/B, M2400A,M2406A, Sonos 2000/2500,tems,M2425A,M2475B	019 	EA035S5A	AHA,	5LD:A,	SNAP,	TPU,	1KΩ	D9			
			EA035C5A	AHA,	5LD:A,	CLIP,	TPU,	1KΩ	D0			
			EA035S3A	AHA,	3LD:A,	SNAP,	TPU,	1KΩ	D7			
			EA035C3A	AHA,	3LD:A,	CLIP,	TPU,	1KΩ	D8			
			EA035S5I	IEC,	5LD:I,	SNAP,	TPU,	1KΩ	D9			
			EA035C5I	IEC,	5LD:I,	CLIP,	TPU,	1KΩ	D0			
EA035S3I			IEC,	3LD:I,	SNAP,	TPU,	1KΩ	D7				
EA035C3I			IEC,	3LD:I,	CLIP,	TPU,	1KΩ	D8				
<b>Primedica</b>	 <b>NEW 2010</b>		EA077S3H	AHA,	3LD:H,	SNAP,	TPU,	10KΩ	I4			
			EA077C3H	AHA,	3LD:H,	CLIP,	TPU,	10KΩ	I4			
			EA077S3IH	IEC,	3LD:IH,	SNAP,	TPU,	10KΩ	I4			
			EA077C3IH	IEC,	3LD:IH,	CLIP,	TPU,	10KΩ	I4			
<b>MEDTRONIC PHYSIO-CONTROL</b>	Lifepak 12 , 120, Lifepak 20		EA066S5A	AHA,	5LD:A,	SNAP,	TPU,	10KΩ	D9			
			EA066C5A	AHA,	5LD:A,	CLIP,	TPU,	10KΩ	D0			
			EA066S3A	AHA,	3LD:H,	SNAP,	TPU,	10KΩ	D7			
			EA066C3A	AHA,	3LD:H,	CLIP,	TPU,	10KΩ	D8			
			EA066S5I	IEC,	5LD:I,	SNAP,	TPU,	10KΩ	D9			
			EA066C5I	IEC,	5LD:I,	CLIP,	TPU,	10KΩ	D0			
			EA066S3I	IEC,	3LD:IH,	SNAP,	TPU,	10KΩ	D7			
			EA066C3I	IEC,	3LD:IH,	CLIP,	TPU,	10KΩ	D8			
	PHYSIO-CONTROL:Lifepak 2,3,33, Scope Monitor II; SIEMENS:Porta Pak 500PB, Porta Pak 500PM; SIMONSEN & WEELS(S&W):CB-600 Defib, DMS-600; SURVIVAL INSTRUMENTS:Safeguard I, II UNIRAD:Series C Echo,VIGGO-SPECTRAMED (Gould):Ekoline 20A, 21, Smith-Kline 20A	282 	EA067S3A	AHA,	3LD:H,	SNAP,	TPU,		B9			
			EA067C3A	AHA,	3LD:H,	CLIP,	TPU,		B0			
			EA067S3I	IEC,	3LD:IH,	SNAP,	TPU,		B9			
			EA067C3I	IEC,	3LD:IH,	CLIP,	TPU,		B0			
			<b>S&amp;W (Artema) Diascope</b>	8000 System, Diascope1,2 , Neoscope, DMS 700, defibrillator, Diascope-NT, Athena Systems, Cardio-Aid, DS Traveller, DSNT, DS1,DS2, MCX/LS		EA018S5A	AHA,	5LD:A,	SNAP,	TPU,		D9
						EA018C5A	AHA,	5LD:A,	CLIP,	TPU,		D0
EA018S3A	AHA,	3LD:H,				SNAP,	TPU,		D7			
EA018C3A	AHA,	3LD:H,				CLIP,	TPU,		D8			
EA018S5I	IEC,	5LD:I,				SNAP,	TPU,		D9			
EA018C5I	IEC,	5LD:I,				CLIP,	TPU,		D0			
EA018S3I	IEC,	3LD:IH,				SNAP,	TPU,		D7			
EA018C3I	IEC,	3LD:IH,				CLIP,	TPU,		D8			
<b>SIEMENS</b>	Sirecust Series 400, 500, 600, Serie 700, 900, 1200, Sirecust 341/610/620/630/720/722 Sirecust 730/732/900/960/961/1260 Sirecust 1281/1801/400/-500 ECG Monitor, Polymed Modul, Modulbox Sirem		EA052S5A	AHA,	5LD:A,	SNAP,	TPU,		B9			
			EA052C5A	AHA,	5LD:A,	CLIP,	TPU,		B0			
			EA052S3A	AHA,	3LD:A,	SNAP,	TPU,		B8			
			EA052C3A	AHA,	3LD:A,	CLIP,	TPU,		B9			
			EA052S5I	IEC,	5LD:I,	SNAP,	TPU,		B9			
			EA052C5I	IEC,	5LD:I,	CLIP,	TPU,		B0			
			EA052S3I	IEC,	3LD:I,	SNAP,	TPU,		B8			
			EA052C3I	IEC,	3LD:I,	CLIP,	TPU,		B9			

Compatible Brand	Compatible Models No.	INSTRUMENT CONNECTOR	REF NO.	STANDARD	LEADS	PATIENT END TERMINAL	Material	Defibrillation withstand	PRICE CODE
<b>SPACELABS</b>	90496 Ultraview, 90369 Ultraview, 90367 Ultraview	 	EA009S5A	AHA,	5LD:A,	SNAP,	TPU,	1KΩ	C9
			EA009C5A	AHA,	5LD:A,	CLIP,	TPU,	1KΩ	C0
			EA009S3A	AHA,	3LD:A,	SNAP,	TPU,	1KΩ	C7
			EA009C3A	AHA,	3LD:A,	CLIP,	TPU,	1KΩ	C8
			EA009S5I	IEC,	5LD:I,	SNAP,	TPU,	1KΩ	C9
			EA009C5I	IEC,	5LD:I,	CLIP,	TPU,	1KΩ	C0
			EA009S3I	IEC,	3LD:I,	SNAP,	TPU,	1KΩ	C7
			EA009C3I	IEC,	3LD:I,	CLIP,	TPU,	1KΩ	C8
<b>SPACELABS BURDICK</b>	<b>GE MEDICAL SYSTEMS:</b> Corometrics 503,Corometrics 504,Corometrics 508; <b>SIEMENS ELECTROMEDICAL SYSTEMS:</b> 806180A Nuclear Camera	 	EA036S3A	AHA,	3LD:H,	SNAP,	TPU,		B9
			EA036C3A	AHA,	3LD:H,	CLIP,	TPU,		B0
			EA036S3I	IEC,	3LD:IH,	SNAP,	TPU,		B9
			EA036C3I	IEC,	3LD:IH,	CLIP,	TPU,		B0
			<b>ADVANCED MEDICAL RESEARCH (AMR):</b> Model 3R; <b>AQUITRON MEDICAL:</b> 9101,9216 (Neonate), <b>ATL ULTRASOUND:</b> Ultramark 8,Ultramark 6,Mark VIII,Mark III,Mark II,600,400; <b>BARD:</b> 500,500 GILSON; <b>BECKMAN (Offner):</b> Offner Dynograph (all),VSM-10C,RM-10 BIOCOM; <b>INC:</b> All models; <b>BOURNS LIFE SYSTEMS:</b> 118-1/2/3; <b>DEL MAR:</b> 2000 Series,2820,2920,2920A; <b>GE MEDICAL SYSTEMS:</b> Critikon Exerstress; <b>MEDICAL ELECTRONICS:</b> All Models (except A30, A40); <b>HEALTHWATCH (Cambridge):</b> MC-4,Simpli Scribe, Versa Scribe, VS-1/2/3; <b>INTERSPEC (Irex) (J&amp;J):</b> System II Echo; <b>KDC:</b> All Models (except Infant IM); <b>MARCONI MEDICAL:</b> 464B, 51-690, CA101, CA105, Echoview 80C, ECO-6, RC-101/M; <b>PIONEER (American Optical):</b> Life Rite 50 ; <b>PPG (Datamedix, Becton-Dickinson):</b> 105 w/lead select, ELD-320W, ELD-400, HR-4, IPA-1/2/3/4, IPA-BD21, IS-5LB, MS-25L, S-5BL , Statscope; <b>PPG (E for M):</b> EAI Series,ECG-6 Series,ECI Series,ECV Series,EDI Series,EEP-8,ERI Series,ETD-8 <b>PPG (Honeywell):</b> 5100 w/ lead select ; <b>SPACELABS BURDICK:</b> CS615,CS625A,CSS20 ,CSS23 ,CSS27,CSS7200,EK1,EK2,EK3,EK4,EK5/A/AS,EK8SS,M200,Medic 1,Medic 2,PA-100,PA-104,PA-105,1600,2100; <b>VIGGO-</b>						
<b>WELCH ALLYN</b>	Propaq 100 Series, Propaq 200 Series, Propaq Encore, Propaq CS , Ultra Smart , Atlas	 	EA021S5A	AHA,	5LD:A,	SNAP,	TPU,	1KΩ	B9
			EA021C5A	AHA,	5LD:A,	CLIP,	TPU,	1KΩ	B0
			EA021S3A	AHA,	3LD:A,	SNAP,	TPU,	1KΩ	B7
			EA021C3A	AHA,	3LD:A,	CLIP,	TPU,	1KΩ	B8
			EA021S5I	IEC,	5LD:I,	SNAP,	TPU,	1KΩ	B9
			EA021C5I	IEC,	5LD:I,	CLIP,	TPU,	1KΩ	B0
			EA021S3I	IEC,	3LD:I,	SNAP,	TPU,	1KΩ	B7
			EA021C3I	IEC,	3LD:I,	CLIP,	TPU,	1KΩ	B8






## Features:


- Ideal for both adult and pediatric patients
- "One-LINE" is compatible with all Infinity patient monitors.
- Fast and easy to apply with clearly labeled,color-coded snaps
- More comfortable for patients,
- Flexible, durable and easy to clean
- "One-LINE" cable design with strategic electrode snap placement
- No need for separate leadwires; one piece design
- Cable can be cleaned or disinfected
- Lightweight wire can be used on a large cross section of patients
- Available in AHA and IEC color coding
- Latex-free


## Type Disignation




Example: E-023P3A indicated Use with PHILIPS M1 series patient monitor, 3-leads, AHA, the Patient end termination is pinch.

## Instrument Connector type





## Instrument Connector



## Yoke type and code

























 Din 3-Ld, Code: <b>C</b>	 Din 5-Ld, Code: <b>C</b>	 Colin, 5-Ld, Code: <b>M</b>
 Lifelinc(LL), 3-Ld, Code: <b>B</b>	 Lifelinc(LL), 5-Ld, Code: <b>B</b>	 TRU-link, 5-Ld, Code: <b>L</b>
 ASP, 3-Ld, Code: <b>H</b>	 ASP, 5-Ld, Code: <b>H</b>	 VS, 5-Ld, Code: <b>E</b>
 AAMI(AA), 3-Ld, Code: <b>D</b>	 AAMI(AA), 5-Ld, Code: <b>D</b>	 NK, 5-Ld, Code: <b>N</b>
 EURO, 3-Ld, Code: <b>J</b>	 EURO, 5-Ld, Code: <b>J</b>	 NK, 6-Ld, Code: <b>P</b>
 Conmed, 3-Ld, Code: <b>G</b>	 Conmed, 5-Ld, Code: <b>G</b>	 Multi-Link, 5Ld, Code: <b>Q</b>
 NEC, 3-Ld, Code: <b>K</b>	 NEC, 3-Ld, Code: <b>K</b>	 Lohmeier multi-link: <b>JM</b>
 Multi-link, Code: <b>R</b>	 Fukuda 16p, Code: <b>Y</b>	 Philips M1668A yoke: <b>DG</b>

Din2.0 Pin: **S**; 2.0Jack yoke: **I**; LL1 type yoke: **U**; NK Din2.0 jack: **V**  
 GE SEER 10Ld: **W**; GE MAX-1: **E**; Philips Trim Yoke: **X**;

























## Type Disignation

E □ □ □ □ □ □

























Color coding: **7A**: 7-Ld, AHA(AAMI); **10A**: 10-Ld, AHA; **7I**: 7-Ld, IEC; **10I**: 10-Ld, IEC; **5A**: 5-Ld, AHA(AAMI); **3H**: 4-Ld, AHA; **3A**: 3-Ld, AHA; **5I**: 5-Ld, IEC; **3IH**: 4-Ld, IEC; **3I**: 3-Ld, IEC; **5AI**: 5-Ld, AHA/IEC; **3H/IH**: 4-Ld, AHA/IEC,  
 Patient end termination type: **S**: snap; **C**: Clip (grabber),  
 If the product is a cable, the lattice is useless, but use the "-" to express, as ED008-3A  
 Brand code, as 041 indicated Bionet; if the product is a leads, the lattice is indicated length, as 040 indicated 40 in.  
 Yoke type: A: one-piece; B: Lifelinc, (LL); C: Din style; D: AA style; E: VS style, other sytle or details as above  
 E: Indicated patient minitor ECG cable and lead wires

Compatible Brand	Compatible Models No.	OEM P/N	REF NO.	Description (Name,Length,Plug,Compatible)	INSTRUMENT CONNECTOR	Yoke (Leads side)	PRICE CODE
<b>Din Style Normal</b>	<b>AIR SHIELDS (NARCO Vickers)</b> :HRRM-71-1/2 ,Porta-Fib I, II, III, System 5 (new), System 6 (Neonatal); <b>ATL ULTRASOUND</b> : 800 ,800 Plus, Apogee , <b>CX ATL ULTRASOUND</b> : CX200, CX400, HDI 9/2000, Interspec Cardiology XL, Mark IV, Mark VII, Ultramark 4 ,Ultramark 4 Plus, Ultramark 5 , Ultramark 7, Ultramark 9; <b>BIOCHEM INTERNATIONAL</b> (BCI):3101,3404,6100,6200, 9100,9200,Advisor ,Autocorr <b>COLIN Medical instruments</b> : BP308, BP408, CBM3000 <b>CRITICARE SYSTEMS, INC.</b> (CSI):504US, 602-4USP, Poet 2, Scholar <b>GE medical systems</b> : Corometrics 556, Critikon 8100T, Critikon Dinamap 8700 Series, Critikon Dinamap 9700 Series <b>IVY biomedical systems</b> :101 (before 1989), 402, 700	D-1340	EC020-3A/I	Din 3Ld Patient Cable,9ft AHA/IEC, straight 6pin>Din 3J, no resistance fit for Din style LeadWires			B9
		D-1540	EC020-5A/I	Din 5Ld Patient Cable,9ft AHA/IEC, straight 6pin>Din 5J, no resistance fit for Din style LeadWires			B0
			EC021-3A/I	Din 3Ld Patient Cable,9ft AHA/IEC, straight 6pin>Din 3J, <b>1KΩ</b> resistance fit for Din style LeadWires			C3
			EC021-5A/I	Din 5Ld Patient Cable,9ft AHA/IEC, straight 6pin>Din 5J, <b>1KΩ</b> resistance fit for Din style LeadWires			C5
<b>Lifelinc (LL) Style Normal</b>	<b>LIFE SCIENCE INSTRUMENTS</b> (LSI):LS5 , 261 <b>MEDICAL DATA electronics</b> (MDE):Escort 100, Escort 101, Escort 102 , Escort 200, Escort 300, Prism , 360SL, 450SL, 521A, PortaPak 80 , PortaPak 90 , PortaPak 1000; <b>Medtronic Physio-control</b> : Lifepak 6 , VSM-1ESF (w/resistors); <b>NELLCOR</b> :N-200 C-LOCK, N-250 C-LOCK, N-1000 , NPB-3900, NPB-4000; <b>NIHON KOHDEN</b> :BSM-8301A, BSM-8302A, BSM-8500A , BSM-8502A, OEC-6102A (s/n>20156), OEC-6105A (s/n>20226) , OEC-7102A , OEC-8108A, TEC-5200A (s/n>20421) , TEC-7100A , TEC-7200A , TEC-7300A; <b>PACE-</b> <b>TECH</b> :Minipack 300 , Minipack 3000, Minipack 3100, MiniMax 4000CL, Vitalmax 800 Plus, Vitalmax 2100, Vitalmax 2200, Vitalmax 4000, Vitalmax 4000CL, Vitalmax 4000 Modular, Vitalmax 4100, Vitalmax 4100CL;		EB020-3A/I	LL 3Ld Patient Cable,9ft AHA/IEC, straight 6pin>LL 3LD, no resistance fit for LL style LeadWires			B9
			EB020-5A/I	LL 5Ld Patient Cable,9ft AHA/IEC, straight 6pin>LL 5LD, no resistance fit for LL style LeadWires			B0
			EB021-3A/I	LL 3Ld Patient Cable,9ft AHA/IEC, straight 6pin>LL 3LD, <b>1KΩ</b> resistance fit for LL style LeadWires			C3
			EB021-5A/I	LL 5Ld Patient Cable,9ft AHA/IEC, straight 6pin>LL 5LD, <b>1KΩ</b> resistance , fit for LL style LeadWires			C5
<b>Conmed Style Normal</b>	<b>PHILIPS MEDICAL SYSTEMS</b> : M3921A A1 Compact Monitor , M3923A A1 Compact Monitor, M3925A A1 Compact Monitor , M3927A A3 Compact Monitor , M3929A A3 Compact Monitor <b>PIONEER</b> (American Optical): 26125 , 33167, 5260 , 5330 , Pulsar 4; <b>PPG</b> (Datamedix, Becton-Dickinson):Guardian 1001, Guardian N1001 <b>SPACELABS BURDICK</b> :Medic 3, Medic 4 , Medic 5 , 4045 Defib, 500 Series, 600 Series, 900 Series (w/o respiration) , 90300 Series, 90400 Series, 90600 Series , 90700 Series, TEK-208, TEK-408, TEK-511 , TEK-512 , TEK-514, TEK-521, PC Bedside , PCMS System; <b>TEKTRONIX (Vitatek)</b> :4045 Defib , 208 , 400 Series, 500 Series; <b>WELCH ALLYN</b> (Protocol):Propaq 100 Series, Propaq 200 Series, Propaq Encore, Propaq CS , Ultra Smart, Atlas <b>ZOLL</b> :D900, PD1200, PD1400		EG020-3A/I	Conmed 3Ld Patient Cable, Length: 9ft ; AHA/IEC, straight 6pin>Conmed 3LD, no resistance, fit for Conmed style LeadWires			C5
			EG020-5A/I	Conmed 5Ld Patient Cable, Length: 9ft ; AHA/IEC, straight 6pin>Conmed 5LD, no resistance, fit for Conmed style LeadWires			C7
			EG021-3A/I	Conmed 3Ld Patient Cable, Length: 9ft ; AHA/IEC, straight 6pin>Conmed 3LD, <b>1KΩ</b> resistance, fit for Conmed style LeadWires			C9
			EG021-5A/I	Conmed 5Ld Patient Cable, Length: 9ft ; AHA/IEC, straight 6pin>Conmed 5LD, <b>1KΩ</b> resistance, fit for Conmed style LeadWires			C0




























## ECG Trunk Cable and Leads for Patient monitor

Compatible Brand	Compatible Models No.	OEM P/N	REF NO.	Description (Name,Length,Plug,Compatible)	INSTRUMENT CONNECTOR	Yoke (Leads side)	PRICE CODE
<b>Bionet Korea</b>	BM3		EB041-3A/I	3Ld Patient Cable,9ft AHA/IEC, Round 8pin>LL-3LD, 1KΩ resistance fit for LL style LeadWires	500 	107 	C2
			EB041-5A/I	5Ld Patient Cable,9ft AHA/IEC, Round 8pin>LL-5LD, 1KΩ resistance fit for LL style LeadWires	015 	105 	C3
<b>Biosys Korea</b>	BPM-103		EC040-3A/I	3Ld Patient Cable,9ft AHA/IEC, Readel 6pin>Din 3J, 1KΩ resistance fit for Din style LeadWires	014 	101 	C3
			EC040-5A/I	5Ld Patient Cable,9ft AHA/IEC, Readel 6pin>Din 5J, 1KΩ resistance fit for Din style LeadWires	014 	102 	C4
<b>Colin</b>	BP88 / BP306		EM062-5I	5Ld Patient Cable,9ft IEC, L type round 6pin>6J, no resistance fit for Colin LeadWires	013 	083 	D2
	BP88S		EM061-5I	5Ld Patient Cable,9ft IEC, straight 6P>6J, no resistance fit for Colin LeadWires	025 	083 	E0
<b>Creative科瑞康 Spurt世博达 Biolight宝莱特</b>	Creative: PC9000S,PC9000,UP6000, Biolight: M12		EC046-3A/I	3Ld Patient Cable,9ft AHA/IEC, Readel 6pin>Din 3J, 1KΩ resistance fit for Din style LeadWires	014 	101 	C3
			EC046-5A/I	5Ld Patient Cable,9ft AHA/IEC, Readel 6pin>Din 5J, 1KΩ resistance fit for Din style LeadWires	014 	102 	C4
<b>Datascope</b>	2000,2001,2002,2100,2200,3000, 3000A ,Ballon Pump,M/D3 w/lead select,M/D3A,M/D4,Passport,Passport IL,Passport XG,Passport XL,Passport 2,Passport 2LT,Spectrum System 90, System 90T , System 96, System 97		EC002-3A/I	Din 3Ld Patient Cable,9ft AHA/IEC, straight 6pin>Din 3J, no resistance fit for Din style LeadWires	005 	101 	B9
			EC002-5A/I	Din 5Ld Patient Cable,9ft AHA/IEC, straight 6pin>Din 5J, no resistance fit for Din style LeadWires	005 	102 	B0
<b>Datex</b>	AS/3, CS/3, S/5, Cardiocap 1, Cardiocap 2 , Cardiocap 5 , Light Monitor , Sattliteplus		EB010-3A/I	3Ld Patient Cable,9ft AHA/IEC, Round 10J>LL 3LD, no resistance fit for LL style LeadWires	023 	107 	C8
			EB010-5A/I	5Ld Patient Cable,9ft AHA/IEC, Round 10J>LL 5LD, no resistance fit for LL style LeadWires	023 	105 	C0

























## ECG Trunk Cable and Leads for Patient monitor

Compatible Brand	Compatible Models No.	OEM P/N	REF NO.	Description (Name,Length,Plug,Compatible)	INSTRUMENT CONNECTOR	Yoke (Leads side)	PRICE CODE
<b>Datex-Ohmeda</b>	Datex-Ohmeda Patient monitor		EH010-3A/I	3Ld Patient Cable,9ft AHA/IEC, Round 10J>ASP 3LD, 4.7KΩ resistance fit for ASP style LeadWires			D1
			EH010-5A/I	5Ld Patient Cable,9ft AHA/IEC, Round 10J>ASP 5LD, no resistance fit for ASP style LeadWires			D2
<b>Drager Siemens</b>	Drager UM1/3/3.1, RM1, PM8010, PM8014, Babyguard 8000, Cato, Babylog8000/8000 plus, PB8800,Serie 8000, Cicero, Parameterbox,		EC057-3A/I	SG 3Ld Patient Cable,9ft AHA/IEC, Round 10pin>Din 3J, no resistance fit for Din style LeadWires			D0
	Multimed-Pot system fur serie SC 6000, SC 6002XL, SC 7000, SC 9000, Art. Nr.3368391 (8.2ft) Art. Nr. 5950196 (4.9ft)	3368391	EQ056-5A/I	SC9000XL Multi-link Cable, 8.2ft AHA/IEC, 1KΩ resistance fit for Euro style LeadWires			O0
<b>Siemens</b>	Sirecust Series 400, 500, 600, Serie 700, 900, 1200, Sirecust 341/610/620/630/720/722, Sirecust 730/732/900/960/961/1260, Sirecust 1281/1801/400/-500, ECG Monitor, Polymed Modul, Modulbox Sirem		EB052-3A/I	3Ld Patient Cable,9ft AHA/IEC, Round 10J>LL 3LD, no resistance fit for LL style LeadWires			B8
			EB052-5A/I	5Ld Patient Cable,9ft AHA/IEC, Round 10J>LL 5LD, no resistance fit for LL style LeadWires			B9
			EJ052-3A/I	3Ld Patient Cable,9ft AHA/IEC, Round 10J>Euro 3LD, no resistance fit for Euro style LeadWires			C5
			EJ052-5A/I	5Ld Patient Cable,9ft AHA/IEC, Round 10J>Euro 5LD, no resistance fit for Euro style LeadWires			C6
<b>GE-Medical Marquette</b>	GEMS:Dash PRO 4000, Dash PRO 3000, Dash PRO 2000, Dash 1000, Eagle Systems, MAC-Lab, Solar 9500, Solar 8000. Solar 3000, Tram 100, Tram 200, Tram 300, Tram 400, Tram 500, Tram 600, Tram 800, Unity Systems		EE007-5A	5Ld Patient Cable,9ft AHA, Round 10J>VS 5LD, no resistance fit for VS style LeadWires			C5
			EE007-5I	5Ld Patient Cable,9ft IEC, Round 10J>VS 5LD, no resistance fit for VS style LeadWires			C8
<b>GE-Hellige</b>	Hellige,Servomed 4 , Generation, SCP 851, SCP 852 , SCP 844, SCS 804,,CardioServ , VICOM-sm		EB011-3A/I	3Ld Patient Cable,9ft AHA/IEC, Round 10P>LL 3LD, no resistance fit for LL style LeadWires			C0
			EB011-5A/I	5Ld Patient Cable,9ft AHA/IEC, Round 10P>LL 5LD, no resistance fit for LL style LeadWires			C0



























## ECG Trunk Cable and Leads for Patient monitor

OEM	OEM MODELS NUMBER	OEM P/N	REF NO.	Description (Name,Length,Plug,Compatible)	INSTRUMENT CONNECTOR	Yoke (Leads side)	PRICE CODE
<b>GE-Hellige</b>	Hellige,Servomed 4 , Generation, SCP 851, SCP 852 , SCP 844, SCS 804,,CardioServ , VICOM-sm		EH011-3A/I	3Ld Patient Cable,9ft AHA/IEC, Round 10J>ASP 3LD, 1KΩ resistance fit for ASP style LeadWires			183 D5
			EH011-5A/I	5Ld Patient Cable,9ft AHA/IEC, Round 10J>ASP 5LD, 1KΩ resistance fit for ASP style LeadWires			221 D6
<b>GE-medical Marquette</b>	Pro1000		EH021-3A/I	3Ld Patient Cable,9ft AHA/IEC, Round 6pin>ASP 3LD, 1KΩ resistance fit for ASP style LeadWires			183 D0
			EH021-5A/I	5Ld Patient Cable,9ft AHA/IEC, Round 6pin>ASP 5LD, 1KΩ resistance fit for ASP style LeadWires			221 E2
<b>Goldway</b>	Goldway: UT4000A,UT4000Apro,UT6000A, UT4000B, UT4000C, UT4000F, UT6000F, UT4000Fpro, UT4800		EB002-3A/I	3Ld Patient Cable,9ft AHA/IEC, Round 6p>LL 3LD, no resistance fit for LL style LeadWires			107 B9
			EB002-5A/I	5Ld Patient Cable,9ft AHA/IEC, Round 6p>LL 5LD, no resistance fit for LL style LeadWires			105 B0
<b>Huntleigh Healthcare</b>	UK MONITOR 		EC047-3A/I	Din 3Ld Patient Cable,9ft AHA/IEC, DB9M>Din 3J, no resistance fit for Din style LeadWires			101 B0
<b>Kontron</b>	7000, 7135, 7135B, 7141, 7210, 7250, 7251, K2000, KAAAT, Kolormon, Micromon, Minimon, Supermon , Vitalmon		EB031-3A/I	3Ld Patient Cable,9ft AHA/IEC, Round 12P>LL 3LD, no resistance fit for LL style LeadWires			107 C4
			EB031-5A/I	5Ld Patient Cable,9ft AHA/IEC, Round 12P>LL 5LD, no resistance fit for LL style LeadWires			105 C5
<b>MEK Korea</b>	MP1000,MP600,MP500		EC029-3A/I	Din 3Ld Patient Cable,9ft AHA/IEC, DB9M>Din 3J, no resistance fit for Din style LeadWires			101 B0
			EC029-5A/I	Din 5Ld Patient Cable,9ft AHA/IEC, DB9M>Din 5J, no resistance fit for Din style LeadWires			102 C1
	M160246, M360697		EB033-3A/I	3Ld Patient Cable,9ft AHA/IEC, MD6M>LL 3LD, no resistance fit for LL style LeadWires			107 C2
			EB033-5A/I	5Ld Patient Cable,9ft AHA/IEC, MD6M>LL 5LD, no resistance fit for LL style LeadWires			105 C3


























## ECG Trunk Cable and Leads for Patient monitor

Compatible Brand	Compatible Models No.	OEM P/N	REF NO.	Description (Name,Length,Plug,Compatible)	INSTRUMENT CONNECTOR	Yoke (Leads side)	PRICE CODE
<b>MINDRAY</b>	PM-7000, PM-8000series, PM-9000series, new PM-6000	0010-30-12378	EJ021-3A/I	3Ld Patient Cable,9ft AHA/IEC, Round 6pin>Euro 3LD, 1KΩ resistance fit for Euro style LeadWires			242 C1
		0010-30-12257	EJ021-5A/I	5Ld Patient Cable,9ft AHA/IEC, Round 6pin>Euro 5LD, 1KΩ resistance fit for Euro style LeadWires			243 C2
	PM-6800,2004年以后的PM-6000, with 12-pin green input connector BeneView (T5, 56, T8) , iPM Series	0010-30-42724	EJ051-3A/I	3Ld Patient Cable,9ft AHA/IEC, straight Round 12P>Euro 3LD, 1KΩ resistance /Inductance fit for Euro style LeadWires			242 D0
		0010-30-42723	EJ051-5A/I	5Ld Patient Cable,9ft AHA/IEC, straight Round 12P>Euro 5LD, 1KΩ resistance /Inductance fit for Euro style LeadWires			243 E1
<b>MINDRAY</b>	PM6000, with 12-pin green input connector	0010-30-12379	EJ050-3A/I	3Ld Patient Cable,9ft AHA/IEC, Type L Round 12P>Euro 3LD, 1KΩ resistance fit for Euro style LeadWires			242 C5
		0010-30-12260	EJ050-5A/I	5Ld Patient Cable,9ft AHA/IEC, Type L Round 12P>Euro 5LD, 1KΩ resistance fit for Euro style LeadWires			243 C6
<b>NEC</b>	3M12	47524	EK053-3I	3Ld Patient Cable,9ft AHA/IEC, Round 16p>NEC 3LD, no resistance fit for 47504 LeadWires			096 D5
		47525	EK053-5I	5Ld Patient Cable,9ft AHA/IEC, Round 16p>NEC 5LD, no resistance fit for 47505 LeadWires			099 D0
<b>Nihon Kohden</b>	OEC-6102A (s/n<20156), OEC-6105A (s/n<20226), OEC-6120A/J/K BSM-7105K		EC037-3A/I	3Ld Patient Cable,9ft AHA/IEC, Round 8p>Din 3LD, no resistance fit for Din style LeadWires			101 C0
		JC-005P	EN037-3I	3Ld Patient Cable,9ft IEC, Round 8p>NK 6J, no resistance fit for BR-019 LeadWires			155 D5
	BSM-7102A,OMP-7201C, OMP-7203C, TEC-5200A (s/n<20421),		EC038-3A/I	3Ld Patient Cable,9ft AHA/IEC, Round 11p>Din 3LD, no resistance fit for Din style LeadWires			101 D5
			EN038-3I	3Ld Patient Cable,9ft IEC, Round 11p>NK 6J no resistance fit for BR-019 LeadWires			155 D0

## ECG Trunk Cable and Leads for Patient monitor







Compatible Brand	Compatible Models No.	OEM P/N	REF NO.	Description (Name,Length,Plug,Compatible)	INSTRUMENT CONNECTOR	Yoke (Leads side)	PRICE CODE
<b>Nihon Kohden</b>	BSM-2301, BSM-2353, BSM-5100, Multi-parameter bedside monitors	JC-906P	EP055-6I	6Ld Patient Cable,9ft ,IEC, Rectangle 12p>NK 8J, no resistance fit for BR-903P, BR-906P			E0
	OPV-1500 Lifescope N	L060 JC-103T	EP080-3I	3Ld Patient Cable,9ft ,IEC, Rectangle 12p>NK 8J, no resistance fit for BR-903P, BR-906P			E0
	BSM-8301A/8302A/8500A,BSM-8502A, OEC-6102A(s/n>20156), OEC-6105A(s/n>20226),OEC-7102A/8108A,TEC-5200A(s/n>20421), 7100A/7200A/7300A;	L152 JC-762V TEC-	EN002-5I	5Ld Patient Cable,9ft IEC, Round 6pin>NK 6J no resistance fit for BR-019 LeadWires			C0
<b>PHILIPS</b>	78352C,7834C, 78834C, M1001A/B, M1002A/B, M1722A/B, M1723A/B, M2475B, M3000A, M3001A, M4735A, VM4/V24	M1500A	ED023-3A	3Ld Patient Cable,9ft AHA, Round 12p>AA-3LD 1KΩ resistance fit for AA style LeadWires			C0
		M1510A	ED023-3I	3Ld Patient Cable,9ft IEC, Round 12p>AA-3LD 1KΩ resistance fit for AA style LeadWires			C0
	78352C,7834C, 78834C, M1001A/B, M1002A/B, M1722A/B, M1723A/B, M2475B, M3000A, M3001A, M4735A, VM4/V24	M1520A	ED023-5A	5Ld Patient Cable,9ft AHA, Round 12p>AA-5LD 1KΩ resistance fit for AA style LeadWires			D4
		M1530A	ED023-5I	5Ld Patient Cable,9ft IEC, Round 12p>AA-5LD 1KΩ resistance fit for AA style LeadWires			D4
<b>PHILIPS</b>	Philips Instruments with 8-Pin Connector: 77000AC,77020A/AC,77020AR,77025A,77030A,77035A,77040A,78208A,78210A,78330A,78331A,78332A,78333A,78341A,78342A,78345A,78346A,78351T,78352A,78353A/B,78354A,78532A/B,78534A/B/C,M1722A/B,M1723A/B,M2400A,M2406A,M2475B,M3500B	M1580A	ED035-3A	3Ld Patient Cable,9ft AHA, Round 8p>AA-3LD 1KΩ resistance fit for AA style LeadWires			C0
		M1590	ED035-3I	3Ld Patient Cable,9ft IEC, Round 8p>AA-3LD 1KΩ resistance fit for AA style LeadWires			C0
	Philips Instruments with 8-Pin Connector: 78208A,78210A,78330A,78331A,78332A,78333A,78341A,78342A,78345A,78346A,78351T,78352A,78353A,78354A,78532A/B,78534A/B/C,M1722A/B,M3500B	M1600A	ED035-5A	5Ld Patient Cable,9ft AHA, Round 8p>AA-5LD 1KΩ resistance fit for AA style LeadWires			D4
		M1610A	ED035-5I	5Ld Patient Cable,9ft IEC, Round 8p>AA-5LD 1KΩ resistance fit for AA style LeadWires			D4
<b>PHILIPS</b>	M3000A, M3001A, M1001A/B,M1002A/B, 78352C,78354C	M1669A	EDG073-3AI	3Ld Patient Cable,9ft IEC/AAMI, Round 12p>3-LD 1KΩ resistance fit for Philips LeadWires			F0
<b>PHILIPS</b>	M3000A, M3001A, M1001A/B,M1002A/B, 78352C,78354C	M1668A	EDG073-5AI	5Ld Patient Cable,9ft IEC/AAMI, Round 12p>5-LD 1KΩ resistance fit for Philips LeadWires			F0









OEM	OEM MODELS NUMBER	OEM P/N	REF NO.	Description (Name,Length,Plug,Compatible)	INSTRUMENT CONNECTOR	Yoke (Leads side)	PRICE CODE	
<b>PHILIPS</b>	78834C 78850A	M1620A	ED036-3A	3Ld Patient Cable,9ft AHA, Round 5p>AA-3LD 1KΩ resistance fit for AA style LeadWires		 086	C0	
			ED036-3I	3Ld Patient Cable,9ft IEC, Round 5p>AA-3LD 1KΩ resistance fit for AA style LeadWires		 086	C0	
<b>S&amp;W (Artema)</b>	8000 System, Diascope1,2 , Neoscope, DMS 700, defibrillator, Diascope-NT, Athena Systems, Cardio-Aid, DS Traveller, DSNT, DS1,DS2, MCX/LS		EB018-3A/I	3Ld Patient Cable,9ft AHA/IEC, Round 10J>LL 3LD, no resistance fit for LL style LeadWires		 107	C8	
			EB018-5A/I	5Ld Patient Cable,9ft AHA/IEC, Round 10J>LL 5LD, no resistance fit for LL style LeadWires		 105	C0	
<b>SPACELABS</b>	90496 Ultraview, 90369 Ultraview,90367 Ultraview	IEC: 700- 0008-07 AHA: 700- 0008-06	EL009-5A/I	5Ld Patient Cable,9ft AHA/IEC, Round 17p>TRU-LINK 5Ld, 1KΩ resistance fit for TRU-LINK LeadWires		 091	C9	
<b>SPACELABS BURDICK</b> Viggo-Spectramed (Statham) WESTINGHOUSE	GE MEDICAL:Corometrics 503,Corometrics 504,Corometrics 508; SIEMENS:806180A Nuclear Camera		EJ036-3A/I	3Ld Patient Cable,9ft AHA/IEC, Round 5pin>Euro 3LD, no resistance fit for Euro style LeadWires		 242	C9	
		<b>ADVANCED MEDICAL RESEARCH (AMR):</b> Model 3R; <b>AEQUITRON MEDICAL:</b> 9101,9216 (Neonate), <b>ATL ULTRASOUND:</b> Ultramark 8,Ultramark 6,Mark VIII,Mark III,Mark II,600,400; <b>BARD:</b> 500,500 GILSON; <b>BECKMAN (Offner):</b> Offner Dynograph (all),VSM-10C,RM-10 BIOCOM; <b>INC:</b> All models; <b>BOURNS LIFE SYSTEMS:</b> 118-1/2/3; <b>DEL MAR:</b> 2000 Series,2820,2920,2920A; <b>GE MEDICAL SYSTEMS:</b> Critikon Exerstress; <b>MEDICAL ELECTRONICS:</b> All Models (except A30, A40); <b>HEALTHWATCH (Cambridge):</b> MC-4,Simpli Scribe, Versa Scribe, VS-1/2/3; <b>INTERSPEC (Irex) (J&amp;J):</b> System II Echo; <b>KDC:</b> All Models (except Infant IM); <b>MARCONI MEDICAL:</b> 464B, 51-690, CA101, CA105, Echoview 80C, ECO-6, RC-101/M; <b>PIONEER (American Optical):</b> Life Rite 50 ; <b>PPG (Datamedix, Becton-Dickinson):</b> 105 w/lead select, ELD-320W, ELD-400, HR-4, IPA-1/2/3/4, IPA-BD21, IS-5LB, MS-25L, S-5BL , Statscope; <b>PPG (E for M):</b> EAI Series,ECG-6 Series,ECI Series,ECV Series,EDI Series,EEP-8,ERI Series,ETD-8 <b>PPG (Honeywell):</b> 5100 w/ lead select ; <b>SPACELABS BURDICK:</b> CS615,CS625A,CSS20 ,CSS23 ,CSS27,CSS7200,EK1,EK2,EK3,EK4,EK5/A/AS,EK8SS,M200,Medic 1,Medic 2,PA-100,PA-104,PA-105,1600,2100; <b>VIGGO-SPECTRAMED (Statham):</b> 6703,Ekoline 21, <b>WESTINGHOUSE:</b> All Models						
<b>Fukuda Denshi</b>	FUKUDA: 7100 patient monitor		EJ032-3A/I	3Ld Patient Cable,9ft AHA/IEC, Round 12P>Euro 3LD, fit for Euro style LeadWires		 242	F2	
				EJ032-5A/I	5Ld Patient Cable,9ft AHA/IEC, Round 12P>Euro 5LD, fit for Euro style LeadWires		 243	F2
		DATASCOPE:Expert FUKUDA DENSHI:DS5100E, DS5300W, DS5400, DS7101/L, DS7141		EY068-5A/I	5Ld and 3Ld Patient Cable,9ft AHA/IEC, Round 12pin> 16 jack non- resistance fit for EY032S3I LeadWires		 373	F2
<b>Lohmeier</b> 			EC078-3A/I	3Ld Patient Cable,9ft <b>IEC/AAMI</b> , 5pin+5Jack plug>Din-3LD non-resistance fit for Din1.5 style Leads		 101	F0	
<b>Lohmeier</b> 	D-81241, original P/N.: 98.060.120/A	98.060.1 20/A	EJM079-5A/I	5Ld Multi-Cable,9ft <b>IEC/AAMI</b> , 27pin plug>Euro 5Ld ECG+2*IBP+SpO2+2*Temp (discharge inside)			O0	

### LEAD WIRES

Low profile, insert molded snap or pinch connectors offer patient comfort and secure attachment to snap electrodes. The snap termination's ergonomic, teardrop design reduces tangling at the patient end connector. Grabber and mini clip terminations have dual sided metal contacts to provide the connection.




Compatible Brand	Productor Guide	OEM P/N	REF NO.	Description (Name,Length,Leads,)	Color coding / Leads	Patient end termination	PRICE CODE		
<b>Bionet</b> <b>Datex</b> <b>Siemens</b> <b>GE-Hellige</b> <b>Kontron</b> <b>Goldway</b> <b>S&amp;W (Artema)</b>   135	suitable for all LL-Plug System , Yoke type: Lifelinc(LL),,Code: <b>B</b>   EB040C3A   EB040C5A		EB024S3A	LL lead wires, <i>24 in. (61 cm), 3-Ld</i>	AHA (AAMI)	Snap	B2		
			EB024C3A		AHA (AAMI)	Grabber (Clip)	B3		
			EB024S3I		IEC	Snap	B2		
			EB024C3I		IEC	Grabber (Clip)	B3		
			EB024S5A	LL lead wires, <i>24 in. (61 cm), 5-Ld</i>	AHA (AAMI)	Snap	B5		
			EB024C5A		AHA (AAMI)	Grabber (Clip)	B6		
			EB024S5I		IEC	Snap	B5		
			EB024C5I		IEC	Grabber (Clip)	B6		
			EB040S3A	LL lead wires, <i>40 in. (102 cm), 3-Ld</i>	AHA (AAMI)	Snap	B2		
			EB040C3A		AHA (AAMI)	Grabber (Clip)	B3		
			EB040S3I		IEC	Snap	B2		
			EB040C3I		IEC	Grabber (Clip)	B3		
			EB040S5A	LL lead wires, <i>40 in. (102 cm), 5-Ld</i>	AHA (AAMI)	Snap	B5		
			EB040C5A		AHA (AAMI)	Grabber (Clip)	B6		
			EB040S5I		IEC	Snap	B5		
			EB040C5I		IEC	Grabber (Clip)	B6		
		<b>Biosys</b> <b>Creative</b> <b>Datascope</b> <b>Nihon Kohden</b> <b>AIR SHIELDS</b> <b>(NARCO, Vickers)</b> <b>ATL</b> <b>ULTRASOUND</b> <b>BIOCHEM</b> <b>INTERNATIONAL</b> <b>COLIN Medical</b> <b>CRITICARE</b> <b>IVY</b>   141	suitable for all Din-Plug System , Yoke type: Din style,Code: <b>C</b>   EC040S3A   EC040S5A		EC018S3A	Din lead wires, <i>18 in. (45 cm), 3-Ld</i>	AHA (AAMI)	Snap	A6
					EC018C3A		AHA (AAMI)	Grabber (Clip)	A7
	EC018S3I			IEC	Snap		A6		
	EC018C3I			IEC	Grabber (Clip)		A7		
	EC018S5A			Din lead wires, <i>18 in. (45 cm), 5-Ld</i>	AHA (AAMI)	Snap	A0		
	EC018C5A				AHA (AAMI)	Grabber (Clip)	B1		
	EC018S5I				IEC	Snap	A0		
	EC018C5I				IEC	Grabber (Clip)	B1		
	EC024S3A			Din lead wires, <i>24 in. (61 cm), 3-Ld</i>	AHA (AAMI)	Snap	A6		
	EC024C3A				AHA (AAMI)	Grabber (Clip)	A7		
	EC024S3I				IEC	Snap	A6		
	EC024C3I				IEC	Grabber (Clip)	A7		
	EC024S5A			Din lead wires, <i>24 in. (61 cm), 5-Ld</i>	AHA (AAMI)	Snap	A0		
	EC024C5A				AHA (AAMI)	Grabber (Clip)	B1		
	EC024S5I				IEC	Snap	A0		
	EC024C5I				IEC	Grabber (Clip)	B1		
	EC040S3A			Din lead wires, <i>40 in. (102 cm), 3-Ld</i>	AHA (AAMI)	Snap	A6		
	EC040C3A				AHA (AAMI)	Grabber (Clip)	A7		
	EC040S3I				IEC	Snap	A6		
	EC040C3I				IEC	Grabber (Clip)	A7		
	EC040S5A	Din lead wires, <i>40 in. (102 cm), 5-Ld</i>	AHA (AAMI)	Snap	A0				
	EC040C5A		AHA (AAMI)	Grabber (Clip)	B1				
	EC040S5I		IEC	Snap	A0				
	EC040C5I		IEC	Grabber (Clip)	B1				


## ECG Trunk Cable and Leads for Patient monitor

Compatible Brand	Productor Guide	OEM P/N	NO.	Description (Name,Length,Leads,)	Color coding / Leads	Patient end termination	PRICE CODE		
<b>PHILIPS</b> M3000A, M3001A, M1001A/B, M1002A/B, 78352C, 78354C	 <b>NEW 2009</b>	M1673A	EDG040S3A	Multi-Link Ldwr, <b>3-Ld</b> Set - Group, <b>40 in</b> (102cm)	AHA (AAMI)	Snap	C0		
		M1674A	EDG040S3I		IEC	Snap	C0		
		M1671A	EDG040C3A		AHA (AAMI)	Grabber (Clip)	C0		
		M1672A	EDG040C3I		IEC	Grabber (Clip)	C0		
		M1644A	EDG040S5A	Multi-Link Ldwr, <b>5-Ld</b> Set - Group, LA, RA, V=40in(102cm) ; LL, RL=61in(155cm) ;	AHA (AAMI)	Snap	D5		
		M1645A	EDG040S5I		IEC	Snap	D5		
		M1968A	EDG040C5A		AHA (AAMI)	Grabber (Clip)	D5		
		M1971A	EDG040C5I		IEC	Grabber (Clip)	D5		
		<b>PHILIPS</b>		M1605A	ED040S3A	AA lead wires, <b>40 in. (102 cm), 3-Ld</b>	AHA (AAMI)	Snap	C5
				M1603A	ED040C3A		AHA (AAMI)	Grabber (Clip)	C6
M1615A	ED040S3I			IEC	Snap		C5		
M1613A	ED040C3I			IEC	Grabber (Clip)		C6		
M1625A	ED040S5A			AA lead wires, <b>5-Ld</b> LA, RA, V 40in(102cm) ; LL, RL 61in(155cm)	AHA (AAMI)	Snap	D4		
M1623A	ED040C5A				AHA (AAMI)	Grabber (Clip)	D5		
M1635A	ED040S5I				IEC	Snap	D4		
M1633A	ED040C5I				IEC	Grabber (Clip)	D5		
 139	suitable for all AA-Plug System , Yoke type: AA style,Code: <b>D</b> <b>OR:For operating room use; includes special filtering; does not allow respiration monitoring.</b>			M1601A	ED040C3A-OR	AA lead wires, <b>OR</b> <b>40 in. (102 cm), 3-Ld</b>	AHA (AAMI)	Grabber (Clip)	C8
				M1611A	ED040C3I-OR	IEC	Grabber (Clip)	C8	
		M1621A	ED040C5A-OR	AA lead wires, <b>5-Ld, OR</b> LA, RA, V 40in(102cm) ; LL, RL 61in(155cm)	AHA (AAMI)	Grabber (Clip)	D6		
		M1631A	ED040C5I-OR	IEC	Grabber (Clip)	D6			
<b>GE--Marquette</b> (Set -Group type Lead wires)	 185	E9008LF	EE029S3A-G	Multi-Link Ldwr Set - Group,, <b>29 in. (76 cm), 3-Ld</b>	AHA (AAMI)	Snap	C1		
		E9003CL	EE029C3A-G		AHA (AAMI)	Grabber (Clip)	C2		
		E9008LH	EE029S3I-G		IEC	Snap	C1		
		E9003CN	EE029C3I-G		IEC	Grabber (Clip)	C2		
		E9008LG	EE051S3A-G	Multi-Link Ldwr Set - Group,, <b>51 in. (130 cm), 3-Ld</b>	AHA (AAMI)	Snap	C1		
		E9003CM	EE051C3A-G		AHA (AAMI)	Grabber (Clip)	C2		
		E9008LJ	EE051S3I-G		IEC	Snap	C1		
		E9003CP	EE051C3I-G		IEC	Grabber (Clip)	C2		
		<b>GE--Marquette</b> (Set -Group type Lead wires)		E9008L	EE029S5A-G	Multi-Link Ldwr Set - Group,, <b>29 in. (76 cm), 5-Ld</b>	AHA (AAMI)	Snap	C1
				E9003CE	EE029C5A-G		AHA (AAMI)	Grabber (Clip)	C2
E9008LB	EE029S5I-G			IEC	Snap		C1		
E9003CG	EE029C5I-G			IEC	Grabber (Clip)		C2		
E9008LA	EE051S5A-G			Multi-Link Ldwr Set - Group,, <b>51 in. (130 cm), 5-Ld</b>	AHA (AAMI)	Snap	C1		
E9003CF	EE051C5A-G				AHA (AAMI)	Grabber (Clip)	C2		
E9008LC	EE051S5I-G				IEC	Snap	C1		
	EE051C5I-G				IEC	Grabber (Clip)	C2		
E9008LD	EE051S5A-GM			Multi-Link Ldwr SetGroup, <b>5-Ld , Mixed Length,</b> RL,LL,V-51in.(130cm) ; RA,LA-29 in.(74cm)	AHA (AAMI)	Snap	C1		
E9003CJ	EE051C5A-GM				AHA (AAMI)	Grabber (Clip)	C2		
E9008LE	EE051S5I-GM	Multi-Link Ldwr SetGroup, <b>5-Ld , Mixed Length,</b> N,F,C-51in. (130cm) , R,L-29 in.(74cm)	IEC	Snap	C1				
	EE051C5I-GM		IEC	Grabber (Clip)	C2				
<b>GE-Medical</b> <b>Marquette</b> (Set -Assembly type Lead wires)	 140	E9008KB	EE029S5A	Multi-Link Ldwr Set <b>29 in. (76 cm), 5-Ld</b>	AHA (AAMI)	Snap	C1		
		E9002ZW	EE029C5A		AHA (AAMI)	Grabber (Clip)	C2		
		E9008KD	EE029S5I		IEC	Snap	C1		
		E9002ZZ	EE029C5I		IEC	Grabber (Clip)	C2		
		E9008KC	EE051S5A	Multi-Link Ldwr Set, <b>51 in. (130 cm), 5-Ld</b>	AHA (AAMI)	Snap	C1		
		E9002ZY	EE051C5A		AHA (AAMI)	Grabber (Clip)	C2		
		E9008KE	EE051S5I		IEC	Snap	C1		
		E9003A	EE051C5I		IEC	Grabber (Clip)	C2		
		E9008KF	EE051S5A-M	Multi-Link Ldwr Set, <b>5-Ld , Mixed Length,</b> RL,LL,V-51in.(130cm); RA,LA-29 in.(74cm)	AHA (AAMI)	Snap	C1		
		E9003AA	EE051C5A-M		AHA (AAMI)	Grabber (Clip)	C2		
	EE051S5I-M	Multi-Link Ldwr Set, <b>5-Ld , Mixed Length,</b> N,F,C-51in.(130cm) ; R,L-29 in.(74cm)	IEC	Snap	C1				
E9003AB	EE051C5I-M		IEC	Grabber (Clip)	C2				










Compatible Brand	Productor Guide	OEM P/N	NO.	Description (Name,Length,Leads,)	Color coding / Leads	Patient end termination	PRICE CODE
<b>Datex PHILIPS PIONEER PPG SPACELABS BURDICK TEKTRONIX (Vitatek) WELCH ALLYN ZOLL</b>   138	suitable for all Conmed-Plug System , Yoke type: Conmed,Code: <u>G</u>  		EG024S3A	Conmed lead wires, <i>24 in. (61 cm), 3-Ld</i>	AHA (AAMI)	Snap	B3
			EG024C3A		AHA (AAMI)	Grabber (Clip)	B4
			EG024S3I		IEC	Snap	B3
			EG024C3I		IEC	Grabber (Clip)	B4
			EG024S5A	Conmed lead wires, <i>24 in. (61 cm), 5-Ld</i>	AHA (AAMI)	Snap	C2
			EG024C5A		AHA (AAMI)	Grabber (Clip)	C3
			EG024S5I		IEC	Snap	C2
			EG024C5I		IEC	Grabber (Clip)	C3
			EG040S3A	Conmed lead wires, <i>40 in. (102 cm), 3-Ld</i>	AHA (AAMI)	Snap	B3
			EG040C3A		AHA (AAMI)	Grabber (Clip)	B4
			EG040S3I		IEC	Snap	B3
			EG040C3I		IEC	Grabber (Clip)	B4
			EG040S5A	Conmed lead wires, <i>40 in. (102 cm), 5-Ld</i>	AHA (AAMI)	Snap	C2
			EG040C5A		AHA (AAMI)	Grabber (Clip)	C3
			EG040S5I		IEC	Snap	C2
			EG040C5I		IEC	Grabber (Clip)	C3
<b>Datex-Ohmeda GE-medical Pro1000</b>   137	suitable for all ASP-Plug System , Yoke type: ASP, Code: <u>H</u>  		EH024S3A	ASP lead wires, <i>24 in. (61 cm), 3-Ld</i>	AHA (AAMI)	Snap	C4
			EH024C3A		AHA (AAMI)	Grabber (Clip)	C5
			EH024S3I		IEC	Snap	C4
			EH024C3I		IEC	Grabber (Clip)	C5
			EH024S5A	ASP lead wires, <i>24 in. (61 cm), 5-Ld</i>	AHA (AAMI)	Snap	D9
			EH024C5A		AHA (AAMI)	Grabber (Clip)	D0
			EH024S5I		IEC	Snap	D9
			EH024C5I		IEC	Grabber (Clip)	D0
			EH040S3A	ASP lead wires, <i>40 in. (102 cm), 3-Ld</i>	AHA (AAMI)	Snap	C4
			EH040C3A		AHA (AAMI)	Grabber (Clip)	C5
			EH040S3I		IEC	Snap	C4
			EH040C3I		IEC	Grabber (Clip)	C5
			EH040S5A	ASP lead wires, <i>40 in. (102 cm), 5-Ld</i>	AHA (AAMI)	Snap	D9
			EH040C5A		AHA (AAMI)	Grabber (Clip)	D0
			EH040S5I		IEC	Snap	D9
			EH040C5I		IEC	Grabber (Clip)	D0
<b>type LM Lead wires</b>   494  			EZ040S3A	LM lead wires, <i>40 in. (102 cm), 3-Ld</i>	AHA (AAMI)	Snap	C4
			EZ040C3A		AHA (AAMI)	Grabber (Clip)	C5
			EZ040S3I		IEC	Snap	C4
			EZ040C3I		IEC	Grabber (Clip)	C5
			EZ040S5A	LM lead wires, <i>40 in. (102 cm), 5-Ld</i>	AHA (AAMI)	Snap	D9
			EZ040C5A		AHA (AAMI)	Grabber (Clip)	D0
			EZ040S5I		IEC	Snap	D9
			EZ040C5I		IEC	Grabber (Clip)	D0
<b>Mindray Siemens SPACELABS BURDICK Viggo- Spectramed (Statham) WESTINGHOUSE</b>   084	suitable for all Euro-Plug System , Yoke type: Euro style, Code: <u>J</u>  		EJ024S3A	Euro style lead wires, <i>24 in. (61 cm), 3-Ld</i>	AHA (AAMI)	Snap	B8
			EJ024C3A		AHA (AAMI)	Grabber (Clip)	B9
			EJ024S3I		IEC	Snap	B8
			EJ024C3I		IEC	Grabber (Clip)	B9
			EJ024S5A	Euro style lead wires, <i>24 in. (61 cm), 5-Ld</i>	AHA (AAMI)	Snap	B0
			EJ024C5A		AHA (AAMI)	Grabber (Clip)	C1
			EJ024S5I		IEC	Snap	B0
			EJ024C5I		IEC	Grabber (Clip)	C1
			EJ040S3A	Euro style lead wires, <i>40 in. (102 cm), 3-Ld</i>	AHA (AAMI)	Snap	B8
			EJ040C3A		AHA (AAMI)	Grabber (Clip)	B9
			EJ040S3I		IEC	Snap	B8
			EJ040C3I		IEC	Grabber (Clip)	B9
			EJ040S5A	Euro style lead wires, <i>40 in. (102 cm), 5-Ld</i>	AHA (AAMI)	Snap	B0
			EJ040C5A		AHA (AAMI)	Grabber (Clip)	C1
			EJ040S5I		IEC	Snap	B0
			EJ040C5I		IEC	Grabber (Clip)	C1

## ECG Trunk Cable and Leads for Patient monitor

Compatible Brand	Productor Guide	OEM P/N	REF NO.	Description (Name,Length,Leads,)	Color coding / Leads	Patient end termination	PRICE CODE
<b>NEC</b> suitable for all NEC-Plug System , Yoke type: NEC style, Code: K	 097		EK029S3I	47504 style lead wires, <u>29 in. (76 cm), 3-Ld</u>	IEC	Snap	B1
		47504	EK029C3I		IEC	Grabber (Clip)	B2
	 100		EK029S5I	47505 style lead wires, <u>29 in. (76 cm), 5-Ld</u>	IEC	Snap	B5
		47505	EK029C5I		IEC	Grabber (Clip)	B6
<b>SPACELABS</b>  166	suitable for all TRU-LINK Plug System , Yoke type: TRU-LINK style, Code: <u>L</u>	700-0007-00	EL024S3A	TUR-LINK style lead wires, <u>24 in. (61 cm), 3-Ld</u>	AHA (AAMI)	Snap	B5
		700-0006-00	EL024C3A		AHA (AAMI)	Grabber (Clip)	B6
		700-0007-01	EL024S3I		IEC	Snap	B5
		700-0006-01	EL024C3I		IEC	Grabber (Clip)	B6
		700-0007-04	EL024S4A	TUR-LINK style lead wires, <u>24 in. (61 cm), 4-Ld</u>	AHA (AAMI)	Snap	B0
		700-0006-04	EL024C4A		AHA (AAMI)	Grabber (Clip)	C1
		700-0007-05	EL024S4I		IEC	Snap	B0
		700-0006-05	EL024C4I		IEC	Grabber (Clip)	C1
		700-0007-08	EL024S5A	TUR-LINK style lead wires, <u>24 in. (61 cm), 5-Ld</u>	AHA (AAMI)	Snap	C4
		700-0006-08	EL024C5A		AHA (AAMI)	Grabber (Clip)	C5
		700-0007-09	EL024S5I		IEC	Snap	C4
		700-0006-09	EL024C5I		IEC	Grabber (Clip)	C5
		700-0007-12	EL024S6A	TUR-LINK style lead wires, <u>24 in. (61 cm), 6-Ld</u>	AHA (AAMI)	Snap	C9
		700-0006-12	EL024C6A		AHA (AAMI)	Grabber (Clip)	C0
		700-0007-13	EL024S6I		IEC	Snap	C9
		700-0006-13	EL024C6I		IEC	Grabber (Clip)	C0
		700-0007-16	EL024S10A	TUR-LINK style lead wires, <u>24 in. (61 cm), 10-Ld</u>	AHA (AAMI)	Snap	E5
		700-0006-16	EL024C10A		AHA (AAMI)	Grabber (Clip)	E6
		700-0007-17	EL024S10I		IEC	Snap	E5
		700-0006-17	EL024C10I		IEC	Grabber (Clip)	E6
		700-0007-02	EL040S3A	TUR-LINK style lead wires, <u>40 in. (102 cm), 3-Ld</u>	AHA (AAMI)	Snap	B5
		700-0006-02	EL040C3A		AHA (AAMI)	Grabber (Clip)	B6
		700-0007-03	EL040S3I		IEC	Snap	B5
		700-0006-03	EL040C3I		IEC	Grabber (Clip)	B6
		700-0007-06	EL040S4A	TUR-LINK style lead wires, <u>40 in. (102 cm), 4-Ld</u>	AHA (AAMI)	Snap	B0
		700-0006-06	EL040C4A		AHA (AAMI)	Grabber (Clip)	C1
		700-0007-07	EL040S4I		IEC	Snap	B0
		700-0006-07	EL040C4I		IEC	Grabber (Clip)	C1
		700-0007-10	EL040S5A	TUR-LINK style lead wires, <u>40 in. (102 cm), 5-Ld</u>	AHA (AAMI)	Snap	C4
		700-0006-10	EL040C5A		AHA (AAMI)	Grabber (Clip)	C5
		700-0007-11	EL040S5I		IEC	Snap	C4
		700-0006-11	EL040C5I		IEC	Grabber (Clip)	C5
		700-0007-14	EL040S6A	TUR-LINK style lead wires, <u>40 in. (102 cm), 6-Ld</u>	AHA (AAMI)	Snap	C9
		700-0006-14	EL040C6A		AHA (AAMI)	Grabber (Clip)	C0
700-0007-15	EL040S6I	IEC	Snap		C9		
700-0006-15	EL040C6I	IEC	Grabber (Clip)		C0		
700-0007-18	EL040S10A	TUR-LINK style lead wires, <u>40 in. (102 cm), 10-Ld</u>	AHA (AAMI)	Snap	E5		
700-0006-18	EL040C10A		AHA (AAMI)	Grabber (Clip)	E6		
700-0007-19	EL040S10I		IEC	Snap	E5		
700-0006-19	EL040C10I		IEC	Grabber (Clip)	E6		



## ECG Trunk Cable and Leads for Patient monitor

Compatible Brand	Productor Guide	OEM P/N	NO.	Description (Name,Length,Leads,)	Color coding / Leads	Patient end termination	PRICE CODE
<b>Colin</b> suitable for all Colin-Plug System , Yoke type: Colin style, Code: <b>M</b>	 009	CL-3601	EM040S3I	COLIN style lead wires, <u>40in. (102 cm), 3-Ld</u>	IEC	Snap	C0
			EM040C3I		IEC	Grabber (Clip)	C0
	 010		EM040S5I	COLIN style lead wires, <u>40in. (102 cm), 5-Ld</u>	IEC	Snap	D5
			EM040C5I		IEC	Grabber (Clip)	D5
<b>Nihon Konden</b> suitable for all NK-5LD Plug System , Yoke type: NK-5LD style, Code: <b>N</b>	 126	BR-019	EN032S3I	NK style lead wires, <u>32 in. 80 cm, 3-Ld</u>	IEC	Snap	C3
			EN032C3I		IEC	Grabber (Clip)	C4
	 125	BR-021P	EN032S5I	NK style lead wires, <u>32 in. 80 cm, 5-Ld</u>	IEC	Snap	C4
			EN032C5I		IEC	Grabber (Clip)	C5
<b>Nihon Konden</b> suitable for all NK-6LD Plug System , Yoke type: NK-6LD style, Code: <b>P</b>	 171/433	BR-913P	EP035S3I	NK style lead wires, <u>35 in. 90 cm, 3-Ld</u>	IEC	Snap	C0
			EP035C3I		IEC	Grabber (Clip)	C0
			EP035S6I	NK style lead wires, <u>35 in. 90 cm, 6-Ld</u>	IEC	Snap	D0
			EP035C6I		IEC	Grabber (Clip)	D0
<b>Nihon Konden</b> suitable for all NK-5LD Plug System , Yoke type: NK-5LD style, Code: <b>R</b>	 124		HA004S3I	NK style lead wires, <u>32 in. 80 cm, 3-Ld</u>	IEC	Snap	B0
		BR-546S	HA004C3I		IEC	Grabber (Clip)	B0
	 123		EV032S3I	NK style lead wires, <u>32 in. 80 cm, 3-Ld</u>	IEC	Snap	B0
		BR-002P	EV032C3I		IEC	Grabber (Clip)	B0
	 122		EV032S5I	NK style lead wires, <u>32 in. 80 cm, 5-Ld</u>	IEC	Snap	C1
		BR-004P	EV032C5I		IEC	Grabber (Clip)	C1
<b>FUKUDA DENSHI</b>	 372		EY032S3I	one-piece lead wires, <u>32 in. 80 cm, 3-Ld</u>	IEC	Snap	C0
			EY032S3A		AHA (AAMI)	Snap	C0

